

Lab for Life.™

Ultra-Waxer 2 | Precision Waxing



Ultra-Waxer 2



NEW: Ultra-Waxer 2

The ultimate electronic waxer just got better, with dual spatulas for increased productivity, tip holder for convenience and organization, and easy to use functions, such as programmable presets, and quick heat function.

Features

- **Systematic.** A complete systematic approach to Dental Prosthetic Waxing.
- **Accurate.** Precise temperature control that is crucial to effective and accurate Waxing.
- **Efficient.** Increased productivity and efficiency are the results of systematic waxing.
- **Proven history.** 30 years as the benchmark of waxing equipment.
- **Comfortable.** Designed with the user's comfort and productivity first and foremost.

Ultra-Spatula and Tips

The Ultra-Spatula features an all aluminum light weight body, with interchangeable Heating Tips, that are color coded for quick identification, foam gripped for ultimate comfort, and ability to change the tip while still hot.



The lost wax casting process begins with a systematic approach to waxing, to ultimately deliver high quality castings every time.

Order Information

34844 **Ultra-Waxer 2**
Contains: Ultra-Tips # 35167, 35168
Ultra-Spatulas # 35177 (2)
Power Cords (US and EU)
Multi Language Instructions

Ultra-Waxer Tips

- 35167 **Small PKT, Blue**
- 35168 **Large PKT, Teal**
- 35169 **Beavertail, Purple**
- 35170 **Small # 7, Grey**
- 35171 **Large # 7, Green**
- 35172 **Needle, Red**
- 35173 **1/2 Hollenback, Black**
- 35174 **Small Denture, Pink**
- 35175 **Large Denture Spoon, Orange**
- 35176 **Denture Spoon, Yellow**
- 35177 **Ultra- Spatula Cord / Handle**

Note: all units are dual voltage, CE and CSA certified