



Maxcem Elite

Maxcem Elite employs the same unique adhesive used in Kerr's OptiBond family of adhesives. The optimized resin matrix and filler systems improve wettability, resulting in greater adhesion and allowing for immediate light curing and better shelf-life stability.

Maxcem Elite's thixotropic property facilitates the placement of material and easy clean up. A proprietary redox initiator system eliminates the inherent discoloration of BPO/tertiary amine initiator systems for a more aesthetic restorations.



24-hour testing. Data from Dr. Carlos Muñoz, SUNY (Buffalo). Available upon request.



24-hour testing. Internal data. Available upon request.